THE TOP IN QUALITY  
for all areas of healthcare

## TOP SYRINGE PUMP TOP-5530/5530RS

On-site Clinical capabilities

A syringe pump with outstanding capabilities and ease of use in the clinical setting, whether on general ward or when administering an anesthetic.

Specifications

|                            |                                                                                                                                                                                                         |                                                                                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Name                       | TOP Syringe Pump                                                                                                                                                                                        |                                                                                                                                                                               |
| Model                      | TOP-5530, TOP-5530RS                                                                                                                                                                                    |                                                                                                                                                                               |
| Power supply               | AC power supply : 100 ~ 127 V ± 10 % 50/60 Hz (for 100 V area)<br>200 ~ 240 V ± 10 % 50/60 Hz (for 200 V area)                                                                                          |                                                                                                                                                                               |
|                            | DC power supply : 12 V ± 10 % (Type: ST55-4)                                                                                                                                                            |                                                                                                                                                                               |
|                            | Internal battery : DC 3.6 V 150 mAh Ni-MH (Type: BP-55)<br>Battery life is approx. 12hours.<br>(At 5 mL/h set using a new battery is fully charged.)                                                    |                                                                                                                                                                               |
|                            | Alkaline battery : LR6 type, AA size alkaline battery DC6V x 4 pcs.<br>Battery life is approx. 24hours.<br>(At 5mL/h set using new batteries)                                                           |                                                                                                                                                                               |
| Input current              | AC Power supply : 0.1 A (for 100V area), 0.05 A (for 200V area)<br>DC power source : 0.4 A                                                                                                              |                                                                                                                                                                               |
| Applicable Syringe         | TOP, B.Braun, TERUMO, B-D, NIPRO, JMS and MONOJECT syringes,<br>Plus one additional syringe registered by user per syringe size<br>10 mL, 20 mL, 30 mL and 50 mL,<br>1% Diprovan Injection-kit 20/50 mL |                                                                                                                                                                               |
| Flow rate range            | mL/h (0.1 mL/h steps)                                                                                                                                                                                   |                                                                                                                                                                               |
|                            | 0.1 ~ 300.0 mL/h                                                                                                                                                                                        | (10 mL syringe)                                                                                                                                                               |
|                            | 0.1 ~ 400.0 mL/h                                                                                                                                                                                        | (20 mL syringe)                                                                                                                                                               |
|                            | 0.1 ~ 500.0 mL/h                                                                                                                                                                                        | (30 mL syringe)                                                                                                                                                               |
| Agent range                | 0.1 ~ 1200.0 mL/h                                                                                                                                                                                       |                                                                                                                                                                               |
|                            | (50 mL syringe, 1% Diprovan Injeciton-kit 20/50 mL)                                                                                                                                                     |                                                                                                                                                                               |
|                            |                                                                                                                                                                                                         |                                                                                                                                                                               |
|                            |                                                                                                                                                                                                         |                                                                                                                                                                               |
| Solution range             | 0.1 ~ 1000.0 mL                                                                                                                                                                                         |                                                                                                                                                                               |
| Weight range               | 0.1 ~ 300.0 kg                                                                                                                                                                                          |                                                                                                                                                                               |
| KVO rate range             | 0.01 ~ 99.99 µg/kg/min                                                                                                                                                                                  |                                                                                                                                                                               |
|                            | 0.01 ~ 99.99 mg/kg/h                                                                                                                                                                                    |                                                                                                                                                                               |
|                            | 0.1 ~ 5.0 mL/h (0.1mL/h steps)                                                                                                                                                                          |                                                                                                                                                                               |
|                            |                                                                                                                                                                                                         |                                                                                                                                                                               |
| Bolus rate                 | 300.0 mL/h (10 mL syringe)                                                                                                                                                                              |                                                                                                                                                                               |
|                            | 400.0 mL/h (20 mL syringe)                                                                                                                                                                              |                                                                                                                                                                               |
|                            | 500.0 mL/h (30 mL syringe)                                                                                                                                                                              |                                                                                                                                                                               |
|                            | 1200.0 mL/h (50 mL syringe, 1% Diprovan Injeciton-kit 20/50 mL)                                                                                                                                         |                                                                                                                                                                               |
| Bolus volume               | mL : 0.1 ~ 50.0 mL (0.1 mL steps)                                                                                                                                                                       |                                                                                                                                                                               |
|                            | µg : 0.1 ~ 50000.0 µg (0.1 µg steps)                                                                                                                                                                    |                                                                                                                                                                               |
|                            | mg : 0.1 ~ 5000.0 mg (0.1 mg steps)                                                                                                                                                                     |                                                                                                                                                                               |
|                            |                                                                                                                                                                                                         |                                                                                                                                                                               |
| Purging                    | Approx. 330 mL/h (10 mL syringe)                                                                                                                                                                        |                                                                                                                                                                               |
|                            | Approx. 620 mL/h (20 mL syringe)                                                                                                                                                                        |                                                                                                                                                                               |
|                            | Approx. 740 mL/h (30 mL syringe)                                                                                                                                                                        |                                                                                                                                                                               |
|                            | Approx. 1200 mL/h (50 mL syringe)                                                                                                                                                                       |                                                                                                                                                                               |
| Volume limit setting       | 0.1 ~ 1000.0 mL                                                                                                                                                                                         |                                                                                                                                                                               |
| Total volume display range | 0.0 ~ 1000.0 mL                                                                                                                                                                                         |                                                                                                                                                                               |
| Infusion accuracy          | Mechanical accuracy : ± 1 %                                                                                                                                                                             |                                                                                                                                                                               |
|                            | Accuracy including syringe : ± 3 %                                                                                                                                                                      |                                                                                                                                                                               |
| Occlusion detection        | 4 steps selection                                                                                                                                                                                       |                                                                                                                                                                               |
|                            | High : 93 ± 33 kPa (700 ± 250 mmHg / 0.95 ± 0.34 kgf/cm2)                                                                                                                                               |                                                                                                                                                                               |
|                            | Medium : 67 ± 27 kPa (500 ± 200 mmHg / 0.68 ± 0.27 kgf/cm2)                                                                                                                                             |                                                                                                                                                                               |
|                            | Low : 40 ± 20 kPa (300 ± 150 mmHg / 0.41 ± 0.20 kgf/cm2)                                                                                                                                                |                                                                                                                                                                               |
|                            | Very low : 20 ± 10 kPa (150 ± 75 mmHg / 0.20 ± 0.10 kgf/cm2)                                                                                                                                            |                                                                                                                                                                               |
|                            | When using TOP Syringe (50 mL) for TOP syringe pump                                                                                                                                                     |                                                                                                                                                                               |
| Overload detection         | Approx. 20N (2.0 kgf) or more (10 mL syringe)                                                                                                                                                           |                                                                                                                                                                               |
|                            | Approx. 40N (4.0 kgf) or more (20 mL syringe)                                                                                                                                                           |                                                                                                                                                                               |
|                            | Approx. 55N (5.5 kgf) or more (30 mL syringe)                                                                                                                                                           |                                                                                                                                                                               |
|                            | Approx. 80N (8.0 kgf) or more (50 mL syringe)                                                                                                                                                           |                                                                                                                                                                               |
| Alarm / Warning            | High priority alarm                                                                                                                                                                                     | Occlusion, Overload, Infusion completion, Battery OFF, Power OFF, Power loss, Malfunction, Syringe barrel, Plunger clamp, Barrel flange, Volume limit completion              |
|                            | Low priority alarm                                                                                                                                                                                      | Low volume, No battery, Low battery, Low alkaline battery, Operation reminder, No rate, No volume limit, No setting value, Sett value check, Prior to volume limit completion |
|                            | Information signal (Warning)                                                                                                                                                                            | Sensor check, Maintenance timer, Standby mode, Power supply switch, Hard limit, Soft limit, Bolus cancel                                                                      |
| Special features           | Power supply switched                                                                                                                                                                                   | The system has entered the build in battery operation mode when power is no longer supplied from the AC/DC power source.                                                      |
|                            | History                                                                                                                                                                                                 | The system stores up to 600/items of operation or operating-state data.                                                                                                       |
|                            | ON/OFF selection                                                                                                                                                                                        | Limits selectable syringe brands or display/hide dose unit, volume limit                                                                                                      |

|                                      |                                                                                |                                                                                                                                                                                                              |
|--------------------------------------|--------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Functions                            | Auto Power OFF                                                                 | During battery operation, when pump is left for 3minutes with infusion stopped or in alarm state, the buzzer sounds. When the buzzer has been left sounding for 3minutes, the power shuts off automatically. |
|                                      | Reminder                                                                       | Unless the cause of alarm is collected within 2 minutes after operator presses the Silence key, alarm will sound again to remind the operator of the problem.                                                |
|                                      | Power supply switch alarm                                                      | Alarm will should when AC power supply has stopped and the unit has switched over to battery operation.                                                                                                      |
|                                      | KVO function                                                                   | After the volume limit is reached, infusion is performed at the set KVO rate for thrombus prevention.                                                                                                        |
|                                      | Standby function                                                               | Reminder alarm can be temporarily reset.                                                                                                                                                                     |
|                                      | Buzzer setting                                                                 | Presence/Absence of wait tone and operation tone can be specified.                                                                                                                                           |
|                                      | History function                                                               | History of Infusion start/stop, alarm, change of flow rate while infusion is under way, power on/off ban be checked.                                                                                         |
|                                      | Key lock function                                                              | During infusion, places the system is Key lock state. In Key lock state, all keys except the power, light and Silence keys are disabled.                                                                     |
|                                      | Syringe brand limit function                                                   | Selectable syringe brands can be limited.                                                                                                                                                                    |
|                                      | Volume limit display Switching function                                        | Whether volume limit is displayed can be specified.                                                                                                                                                          |
|                                      | Battery refresh function                                                       | Performs refresh operation in order to maintain the battery performance.                                                                                                                                     |
|                                      | Maintenance timer function                                                     | When the preset timing has been reached, a screen to prompt regular inspection is displayed.                                                                                                                 |
| Operating conditions                 | Flow rate upper limit value setting function                                   | Flow rate upper limit value can be set.                                                                                                                                                                      |
|                                      | User syringe registration function                                             | A syringe not registered to this machine can be registered.                                                                                                                                                  |
| Transportation and Storage condition | Ambient temperature : 5 to 40 °C (While charging the internal battery: 5-30°C) |                                                                                                                                                                                                              |
|                                      | Relative humidity : 20 to 90 % (Noncondensing)                                 |                                                                                                                                                                                                              |
| Service life                         | Atmospheric pressure : 70 to 106 kPa                                           |                                                                                                                                                                                                              |
|                                      | Ambient temperature : -20 to 45 °C                                             |                                                                                                                                                                                                              |
| Classification                       | Relative humidity : 10 to 95 % (Noncondensing)                                 |                                                                                                                                                                                                              |
|                                      | Atmospheric pressure : 50 to 106 kPa                                           |                                                                                                                                                                                                              |
| Fuse                                 | T0.2AL (for 100 V area), T0.1AL (for 200 V area)                               |                                                                                                                                                                                                              |
| Detentions                           | 320(W) x 90(H) x 160(D) mm                                                     |                                                                                                                                                                                                              |
| Weight                               | Approx. 2.0 kg                                                                 |                                                                                                                                                                                                              |
| Accessories                          | AC power cable x1 Operation guide x1, Administrator Manual x1                  |                                                                                                                                                                                                              |
| Material used                        | Case                                                                           | ABS (Acrylonitrile Butadiene Styrene)                                                                                                                                                                        |
|                                      | Syringe barrel clamp                                                           | Glass fiber-reinforced Polycarbonate (30%)                                                                                                                                                                   |
| Options                              | Syringe plunger slider                                                         | Glass fiber-reinforced Polycarbonate (30%)                                                                                                                                                                   |
|                                      | Operation panel                                                                | Polyester film                                                                                                                                                                                               |
|                                      | LED cover                                                                      | Methyl methacrylate polymer                                                                                                                                                                                  |
|                                      | Pole clamp, Multiple pump mount, Infusino stand                                |                                                                                                                                                                                                              |

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**MEDITOP CORPORATION (MALAYSIA) SDN. BHD.**

yang beralamat di:  
which is located at:

**NO. 3, PERSIARAN USAHAWAN, TAMAN IKS,  
SEKSYEN 9, ,  
43650  
BANDAR BARU BANGI SELANGOR DARUL  
EHSAN**

bagi mengesahkan peranti perubatan seperti yang dinyatakan dalam Lampiran 1 adalah berdaftar di bawah Seksyen 5(1) Akta 737.  
to confirm that the medical device as detailed out in Attachment 1 is registered under Section 5(1) of Act 737.

Pendaftaran ini diberikan tertakluk kepada peruntukan-peruntukan di bawah Akta 737 dan peraturan-peraturan yang dibuat dibawahnya serta syarat-syarat seperti di Lampiran 2.  
This registration is granted subject to the provisions under Act 737 and its subsidiary legislations and the conditions as in Attachment 2.

**AHMAD SHARIFF BIN HAMBALI**  
KETUA EKSEKUTIF  
CHIEF EXECUTIVE  
**PIHAK BERKUASA PERANTI PERUBATAN**  
MEDICAL DEVICE AUTHORITY

No. Pendaftaran: **GC4581638216**  
Registration No.:

Butir-butir peranti perubatan yang didaftarkan  
Particulars of the registered medical device

Nama Peranti Perubatan **SYRINGE PUMP**  
Medical Device Name

Kelas **CLASS C** Jenama **TOP**  
Class Brand

Kelompok **FAMILY**  
Group

Nama dan alamat pembuat: **TOP CORPORATION**  
Name and address of manufacturer **19-10 SENJU NAKAI-CHO, ADACHI-KU, TOKYO 120-0035, JAPAN, JAPAN**

APPENDIX

| NO | NAME AS PER DEVICE LABEL | IDENTIFIER | BRIEF DESCRIPTION OF ITEM                                                                                                                                                                                                                                                               |
|----|--------------------------|------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1  | TOP SYRINGE PUMP         | TOP-5520RS | These are a micro-continuous infusion pump for parenteral nutrition and blood transfusion to immature infants and neonates, infusion of chemotherapeutic agents such as antineoplastic, oxytocic, anticoagulant, and anesthetic agents in the Intensive Care Unit / Critical Care Unit. |
| 2  | TOP SYRINGE PUMP         | TOP-5530   | These are a micro-continuous infusion pump for parenteral nutrition and blood transfusion to immature infants and neonates, infusion of chemotherapeutic agents such as antineoplastic, oxytocic, anticoagulant, and anesthetic agents in the Intensive Care Unit / Critical Care Unit. |

| NO                    | NAME AS PER DEVICE LABEL | IDENTIFIER | BRIEF DESCRIPTION OF ITEM                                                                                                                                                                                                                                                               |
|-----------------------|--------------------------|------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 3                     | TOP SYRINGE PUMP         | TOP-5530RS | These are a micro-continuous infusion pump for parenteral nutrition and blood transfusion to immature infants and neonates, infusion of chemotherapeutic agents such as antineoplastic, oxytocic, anticoagulant, and anesthetic agents in the Intensive Care Unit / Critical Care Unit. |
| "End Of Product List" |                          |            |                                                                                                                                                                                                                                                                                         |

**SYARAT – SYARAT PENDAFTARAN PERANTI PERUBATAN**  
**MEDICAL DEVICE REGISTRATION CONDITIONS**

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## **2.0 PINDAAN PENDAFTARAN PERANTI PERUBATAN AMENDMENT OF MEDICAL DEVICE REGISTRATION**

- 2.1 Sebarang pindaan kepada maklumat yang berkaitan peranti perubatan yang berdaftar hendaklah dimaklumkan kepada Pihak Berkuasa secara rasmi mengikut garis panduan yang ditetapkan oleh Pihak Berkuasa. Pihak Berkuasa berhak memberikan kelulusan atau menolak permohonan pindaan tersebut.  
*Any amendments to the information concerning registered medical device shall be notified to the Authority in accordance to the guidelines set by the Authority. The Authority reserves the right to grant approval or reject the application for such amendments.*

## **3.0 PEMBATALAN SIJIL PENDAFTARAN PERANTI PERUBATAN CANCELLATION OF MEDICAL DEVICE REGISTRATION CERTIFICATE**

- 3.1 Sijil Pendaftaran Peranti Perubatan boleh dibatalkan seperti yang dinyatakan dalam Seksyen 9, Akta 737.  
*Medical Device Registration may be cancelled as prescribed in Section 9 of Act 737.*
- 3.2 Mana-mana peranti perubatan yang dibatalkan Sijil Pendaftarannya, tidak boleh diimport, dieksport atau diletakkan dalam pasaran.  
*Any Medical Device which the registration certificate has been cancelled shall not be imported, exported or placed in the market.*

**4.0 HAK PIHAK BERKUASA**  
**THE AUTHORITY OWNERSHIP**

- 4.1 Sijil Pendaftaran Peranti Perubatan yang dikeluarkan secara manual atau elektronik adalah **Hak Milik Pihak Berkuasa**.  
*The Authority retains the ownership of every Medical Device Registration Certificate issued by any means.*
- 4.2 Sekiranya berlaku kehilangan atau kerosakan Sijil Pendaftaran Peranti Perubatan, hendaklah dimaklumkan kepada Pihak Berkuasa dan setiap penggantian sijil akan dikenakan caj perkhidmatan.  
*Any loss or damage to the Medical Device Registration Certificate shall be notified to the Authority and every replacement of certificate shall be liable with service charge rendered.*

**5.0 TUGAS DAN TANGGUNGJAWAB**  
**ROLES AND RESPONSIBILITIES**

- 5.1 Establismen hendaklah mematuhi Akta 737, peraturan-peraturan di bawah Akta dan syarat-syarat Pendaftaran Peranti Perubatan.  
*Establishment shall comply with Act 737, its subsidiary regulations and registration conditions.*