THE TOP IN QUALITY
for all areas of healthcare

TOP SYRINGE PUMP TOP-5530/5530RS

On-site Clinical capabilities

A syringe pump with outstanding capabilities and ease of use in the clinical setting, whether on general ward or when administering an anesthetic.

Specifications

Name	TOP Syringe Pump	
Model	TOP-5530, TOP-5530RS	
Power supply	AC power supply : 100 ~ 127 V ± 10 % 50/60 Hz (for 100 V area) 200 ~ 240 V ± 10 % 50/60 Hz (for 200 V area)	
	DC power supply : 12 V ± 10 % (Type: ST55-4) Internal battery : DC 3.6 V 150 mAh Ni-MH (Type: BP-55) Battery life is approx. 12hours. (At 5 mL/h set using a new battery is fully charged.)	
	Alkaline battery : LR6 type, AA size alkaline battery DC6V x 4 pcs. Battery life is approx. 24hours. (At 5mL/h set using new batteries)	
Input current	AC Power supply : 0.1 A (for 100V area), 0.05 A (for 200V area) DC power source : 0.4 A	
Applicable Syringe	TOP, B.Braun, TERUMO, B-D, NIPRO, JMS and MONOJECT syringes, Plus one additional syringe registered by user per syringe size 10 mL, 20 mL, 30 mL and 50 mL, 1% Diprovan Injection-kit 20/50 mL	
Flow rate range	mL/h (0.1 mL/h steps)	
	0.1 ~ 300.0 mL/h	(10 mL syringe)
	0.1 ~ 400.0 mL/h	(20 mL syringe)
	0.1 ~ 500.0 mL/h	(30 mL syringe)
Agent range	0.1 ~ 1200.0 mL/h (50 mL syringe, 1% Diprovan Injeciton-kit 20/50 mL)	
Solution range	0.1 ~ 1000.0 mg	
Weight range	0.1 ~ 300.0 kg	
KVO rate range	0.01 ~ 99.99 µg/kg/min	
	0.01 ~ 99.99 mg/kg/h	
	0.1 ~ 5.0 mL/h (0.1mL/h steps)	
Bolus rate	300.0 mL/h	(10 mL syringe)
	400.0 mL/h	(20 mL syringe)
	500.0 mL/h	(30 mL syringe)
	1200.0 mL/h	(50 mL syringe, 1% Diprovan Injeciton-kit 20/50 mL)
Bolus volume	mL : 0.1 ~ 50.0 mL	(0.1 mL steps)
	µg : 0.1 ~ 50000.0 µg	(0.1 µg steps)
	mg : 0.1 ~ 5000.0 mg	(0.1 mg steps)
Purging	Approx. 330 mL/h	(10 mL syringe)
	Approx. 620 mL/h	(20 mL syringe)
	Approx. 740 mL/h	(30 mL syringe)
	Approx. 1200 mL/h	(50 mL syringe)
Volume limit setting	0.1 ~ 1000.0 mL	
Total volume display range	0.0 ~ 1000.0 mL	
Infusion accuracy	Mechanical accuracy : ± 1 %	
	Accuracy including syringe : ± 3 %	
Occlusion detection	4 steps selection	
	High	93 ± 33 kPa (700 ± 250 mmHg / 0.95 ± 0.34 kgf/cm ²)
	Medium	67 ± 27 kPa (500 ± 200 mmHg / 0.68 ± 0.27 kgf/cm ²)
	Low	40 ± 20 kPa (300 ± 150 mmHg / 0.41 ± 0.20 kgf/cm ²)
	Very low	20 ± 10 kPa (150 ± 75 mmHg / 0.20 ± 0.10 kgf/cm ²)
	When using TOP Syringe (50 mL) for TOP syringe pump	
Overload detection	Approx. 20N (2.0 kgf) or more	(10 mL syringe)
	Approx. 40N (4.0 kgf) or more	(20 mL syringe)
	Approx. 55N (5.5 kgf) or more	(30 mL syringe)
	Approx. 80N (8.0 kgf) or more	(50 mL syringe)
Alarm / Warning	High priority alarm	Occlusion, Overload, Infusion completion, Battery OFF, Power OFF, Power loss, Malfunction, Syringe barrel, Plunger clamp, Barrel flange, Volume limit completion
	Low priority alarm	Low volume, No battery, Low battery, Low alkaline battery, Operation reminder, No rate, No volume limit, No setting value, Sett value check, Prior to volume limit completion
	Information signal (Warning)	Sensor check, Maintenance timer, Standby mode, Power supply switch, Hard limit, Soft limit, Bolus cancel
Special features	Power supply switched	The system has entered the build in battery operation mode when power is no longer supplied from the AC/DC power source.
	History	The system stores up to 600items of operation or operating-state data.
	ON/OFF selection	Limits selectable syringe brands or display/hide dose unit, volume limit

Functions	Auto Power OFF	During battery operation, when pump is left for 3minutes with infusion stopped or in alarm state, the buzzer sounds. When the buzzer has been left sounding for 3minutes, the power shuts off automatically.
	Reminder	Unless the cause of alarm is collected within 2 minutes after operator presses the Silence key, alarm will sound again to remind the operator of the problem.
	Power supply switch alarm	Alarm will should when AC power supply has stopped and the unit has switched over to battery operation.
	KVO function	After the volume limit is reached, infusion is performed at the set KVO rate for thrombus prevention.
	Standby function	Reminder alarm can be temporarily reset.
	Buzzer setting	Presence/Absence of wait tone and operation tone can be specified.
	History function	History of Infusion start/stop, alarm, change of flow rate while infusion is under way, power on/off can be checked.
	Key lock function	During infusion, places the system is Key lock state. In Key lock state, all keys except the power, light and Silence keys are disabled.
	Syringe brand limit function	Selectable syringe brands can be limited.
	Volume limit display	Whether volume limit is displayed can be specified.
	Switching function	Performs refresh operation in order to maintain the battery performance.
	Battery refresh function	When the preset timing has been reached, a screen to prompt regular inspection is displayed.
Operating conditions	Flow rate upper limit value setting function	Flow rate upper limit value can be set.
	User syringe registration function	A syringe not registered to this machine can be registered.
Transportation and Storage condition	Ambient temperature : 5 to 40 °C (While charging the internal battery: 5-30°C)	
	Relative humidity : 20 to 90 % (Noncondensing)	
Service life	Atmospheric pressure : 70 to 106 kPa	
	Ambient temperature : -20 to 45 °C	
Classification	Relative humidity : 10 to 95 % (Noncondensing)	
	Atmospheric pressure : 50 to 106 kPa	
Fuse	6years (Based on self certification using TOP corporation data)	
Detentions	Class II  and Internally powered equipment, Type CF  , IPX2 (Drip proof)	
Weight	T0.2AL (for 100 V area), T0.1AL (for 200 V area)	
Accessories	320(W) x 90(H) x 160(D) mm	
Material used	Approx. 2.0 kg	
	AC power cable x1	Operation guide x1, Administrator Manual x1
Options	Case	ABS (Acrylonitrile Butadiene Styrene)
	Syringe barrel clamp	Glass fiber-reinforced Polycarbonate (30%)
	Syringe plunger slider	Glass fiber-reinforced Polycarbonate (30%)
	Operation panel	Polyester film
	LED cover	Methyl methacrylate polymer
	Pole clamp, Multiple pump mount, Infusino stand	

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MEDITOP CORPORATION (MALAYSIA) SDN. BHD.

yang beralamat di:
which is located at:

**NO. 3, PERSIARAN USAHAWAN, TAMAN IKS,
SEKSYEN 9, ,
43650
BANDAR BARU BANGI SELANGOR DARUL
EHSAN**

bagi mengesahkan peranti perubatan seperti yang dinyatakan dalam Lampiran 1 adalah berdaftar di bawah Seksyen 5(1) Akta 737.
to confirm that the medical device as detailed out in Attachment 1 is registered under Section 5(1) of Act 737.

Pendaftaran ini diberikan tertakluk kepada peruntukan-peruntukan di bawah Akta 737 dan peraturan-peraturan yang dibuat dibawahnya serta syarat-syarat seperti di Lampiran 2.
This registration is granted subject to the provisions under Act 737 and its subsidiary legislations and the conditions as in Attachment 2.

AHMAD SHARIFF BIN HAMBALI
KETUA EKSEKUTIF
CHIEF EXECUTIVE
PIHAK BERKUASA PERANTI PERUBATAN
MEDICAL DEVICE AUTHORITY

No. Pendaftaran: **GC4581638216**
Registration No.:

Butir-butir peranti perubatan yang didaftarkan
Particulars of the registered medical device

Nama Peranti Perubatan **SYRINGE PUMP**
Medical Device Name

Kelas **CLASS C** Jenama **TOP**
Class Brand

Kelompok **FAMILY**
Group

Nama dan alamat pembuat: **TOP CORPORATION**
Name and address of manufacturer **19-10 SENJU NAKAI-CHO, ADACHI-KU, TOKYO 120-0035, JAPAN, JAPAN**

APPENDIX

NO	NAME AS PER DEVICE LABEL	IDENTIFIER	BRIEF DESCRIPTION OF ITEM
1	TOP SYRINGE PUMP	TOP-5520RS	These are a micro-continuous infusion pump for parenteral nutrition and blood transfusion to immature infants and neonates, infusion of chemotherapeutic agents such as antineoplastic, oxytocic, anticoagulant, and anesthetic agents in the Intensive Care Unit / Critical Care Unit.
2	TOP SYRINGE PUMP	TOP-5530	These are a micro-continuous infusion pump for parenteral nutrition and blood transfusion to immature infants and neonates, infusion of chemotherapeutic agents such as antineoplastic, oxytocic, anticoagulant, and anesthetic agents in the Intensive Care Unit / Critical Care Unit.

NO	NAME AS PER DEVICE LABEL	IDENTIFIER	BRIEF DESCRIPTION OF ITEM
3	TOP SYRINGE PUMP	TOP-5530RS	These are a micro-continuous infusion pump for parenteral nutrition and blood transfusion to immature infants and neonates, infusion of chemotherapeutic agents such as antineoplastic, oxytocic, anticoagulant, and anesthetic agents in the Intensive Care Unit / Critical Care Unit.
"End Of Product List"			

SYARAT – SYARAT PENDAFTARAN PERANTI PERUBATAN
MEDICAL DEVICE REGISTRATION CONDITIONS

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2.0 PINDAAN PENDAFTARAN PERANTI PERUBATAN AMENDMENT OF MEDICAL DEVICE REGISTRATION

- 2.1 Sebarang pindaan kepada maklumat yang berkaitan peranti perubatan yang berdaftar hendaklah dimaklumkan kepada Pihak Berkuasa secara rasmi mengikut garis panduan yang ditetapkan oleh Pihak Berkuasa. Pihak Berkuasa berhak memberikan kelulusan atau menolak permohonan pindaan tersebut.
Any amendments to the information concerning registered medical device shall be notified to the Authority in accordance to the guidelines set by the Authority. The Authority reserves the right to grant approval or reject the application for such amendments.

3.0 PEMBATALAN SIJIL PENDAFTARAN PERANTI PERUBATAN CANCELLATION OF MEDICAL DEVICE REGISTRATION CERTIFICATE

- 3.1 Sijil Pendaftaran Peranti Perubatan boleh dibatalkan seperti yang dinyatakan dalam Seksyen 9, Akta 737.
Medical Device Registration may be cancelled as prescribed in Section 9 of Act 737.
- 3.2 Mana-mana peranti perubatan yang dibatalkan Sijil Pendaftarannya, tidak boleh diimport, dieksport atau diletakkan dalam pasaran.
Any Medical Device which the registration certificate has been cancelled shall not be imported, exported or placed in the market.

4.0 HAK PIHAK BERKUASA
THE AUTHORITY OWNERSHIP

- 4.1 Sijil Pendaftaran Peranti Perubatan yang dikeluarkan secara manual atau elektronik adalah **Hak Milik Pihak Berkuasa**.
The Authority retains the ownership of every Medical Device Registration Certificate issued by any means.
- 4.2 Sekiranya berlaku kehilangan atau kerosakan Sijil Pendaftaran Peranti Perubatan, hendaklah dimaklumkan kepada Pihak Berkuasa dan setiap penggantian sijil akan dikenakan caj perkhidmatan.
Any loss or damage to the Medical Device Registration Certificate shall be notified to the Authority and every replacement of certificate shall be liable with service charge rendered.

5.0 TUGAS DAN TANGGUNGJAWAB
ROLES AND RESPONSIBILITIES

- 5.1 Establismen hendaklah mematuhi Akta 737, peraturan-peraturan di bawah Akta dan syarat-syarat Pendaftaran Peranti Perubatan.
Establishment shall comply with Act 737, its subsidiary regulations and registration conditions.